

Order form - select the correct Cover

 Flex Neo AK

 XS COVER

Cover circumference in [cm]


Min. 38 for Genium, GeniumX3, C-Leg4, Rheo
Min. 40 for GeniumX4

36
38
40
42

Reduce by 5% to 15%
For example :
41cm - 5% = 38,9cm

38

The height of the knee in the range of 31-55 in [cm]



Foot manufacturer, type and size

For example : OB Taleo 1C50 L29

Knee manufacturer and type

For example : OB 3R80

Adapter

- 30
 34

e.g.:
4WR95=3

Matte and semi-matte colours - tattoos unavailable

- | | | | |
|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| <input type="checkbox"/> 21 | <input type="checkbox"/> 41 | <input type="checkbox"/> 61 | <input type="checkbox"/> 81 |
| <input type="checkbox"/> 23 | <input type="checkbox"/> 43 | <input type="checkbox"/> 63 | <input type="checkbox"/> 83 |
| <input type="checkbox"/> 25 | <input type="checkbox"/> 45 | <input type="checkbox"/> 65 | <input type="checkbox"/> 85 |

Shiny colours & matt varnish

- | | | |
|-------------------------------------|---|-------------------------------------|
| <input type="checkbox"/> Skin Tone1 | <input type="checkbox"/> Skin Tone2 | <input type="checkbox"/> Skin Tone3 |
| <input type="checkbox"/> Graphite | <input type="checkbox"/> Grey | <input type="checkbox"/> White |
| <input type="checkbox"/> Black | <input type="checkbox"/> Matt varnish (charged extra) | |

Tattoo (charged extra)

Available only for Skin Tone3, White and Grey colour


- Catalog number
e.g.: 08-1001

- Own


Check list :

4 Photos The camera must be placed perpendicular to the legs, in a low position, and both legs must be centered on the photo. Please remember that all prosthetic components must be completely visible, including the socket. Please remove shoes, socks and any object that covers the prosthesis.

Patients do not need to remove the socks from the healthy leg.



VERY IMPORTANT



Please place your phone no higher than 20 cm above the floor !!!!!!!!!!!!!!!



- All our covers can be shortened from the bottom by several mm.
- In the factory state, the lower part of the cover is flat.
- In the case of the foot as shown in the picture, the lower part of the cover can be ground to the shape of the foot.

COLOURS - do not send.

Matte and semi-matte colours - tattoos unavailable



Shiny colours - tattoos available

