



# FLEX NEO AK


ABOVE KNEE AMPUTATION, TRANSFEMORAL



## TECHNICAL DETAILS



36; 38; 40; 42 cm



Any in the range of 32-53 cm  
\* covers can be shortened by a dozen mm




- Available for all knee joints.
- Easy and fast to put on thanks to magnets.
- Waterproof.
- Calf circumference: 36; 38; 40; 42 cm.
- Height: any from the 32-53cm range.
- Weight: approx. 500 g.
- Shape: anatomical
- Material: soft and flexible, it is highly resistant to water and bacteria.
- The cross section of the lower part of the cover is elliptical.
- Painless Tattoos - extra paid option.
- Matte version - extra paid option.
- There is a gap of approximately 7 mm between the cover and the foot. For very soft feet, such as carbon fiber feet, the distance should be significantly greater than 7 mm.


## COMPARISON

	Flex Neo AK	Flex Neo MAX AK
<b>MAGNETS</b>	Yes	Yes
<b>HEIGHT</b>	Any of the range of 32 - 53 cm	Any
<b>CIRCUMFERENCE</b>	36; 38; 40; 42 cm	Any
<b>SHAPE</b>	Anatomical (left or right)	Reconstruction of calf based on a scan.
<b>A SCAN IS REQUIRED TO ORDER</b>	No	Yes


# COLOURS




**Skin tone**





**Grey**




**White**



**Black**



**Graphite**



**Glossy**

**Matte**