




# KIDS


## ABOVE & BELOW KNEE AMPUTATION



# TECHNICAL DETAILS




Any





Any

\* covers can be shortened by a dozen mm




Any





The cross-section of the lower part of the cover is similar to the shape of the foot shell.

- Waterproof.
- Calf circumference: any.
- Height: any.
- Weight: approx. 200 g.
- Shape: anatomical.
- Material: soft and flexible, it is highly resistant to water and bacteria.
- Painless Tattoos - extra paid option.
- Matte version - extra paid option.


# COLOURS




Skin tone





Grey




White



Black



Graphite



Glossy

Matte