


Order Form - Flex **Neo** AK Transfemoral


Cover circumference in [cm]

- 36
- 38
- 40
- 42



We recommend reducing the FlexCover circumference by a few percent, and even above ten percent for a very large calf. According to our customers' opinion, the FlexCover looks good with a circumference of less than 40 cm

The height of the knee in the range of 32-53 cm in [cm]



We measure from the level of the foot shell to the top of the knee. The gap between the cover and the foot will be taken into account by us, depending on the type of foot.

Additionally, the covers can be shortened by a few cm. In the case of a parabolic foot shell, the height measurement must take into account the lowest point. In this case, the cover should be ground by yourself to the shape of the foot.



For example : OB 3R80

Knee manufacturer and type

Adapter

- 30
- 34

Foot

- Left
- Right

Colour

- SkinTone
- Black
- Grey
- White


Finish

- Matt varnish
(charged extra)

Tattoo (charged extra)

- Catalog number from
the our website.
e.g.; 08-1001

Only for SkinTone, Red and White colour



- All our covers can be shortened from the bottom by a few cm.
- In the delivered condition, the lower part of the cover is flat.
- In the case of a parabolic foot shell, the lower part of the cover can be ground down to the shape of the foot.
- Inner dimensions of the lower part of the cover: 58 mm X 85 mm.

