


Prosthetic socket circumference at the very end of the cover in [mm]





To take the measurement correctly, please mark the height of the cover in several places and place a ruler below the marking. The measurement should be made very carefully because moving the ruler even by a few millimetres up or down causes an incorrect reading of the circumference even by several centimetres.



Height in [mm]

We measure from the level of the foot shell to the top of the knee. The gap between the cover and the foot will be taken into account by us, depending on the type of foot.



A [mm]

B [mm]

Foot manufacturer, type and size

For example : OB Taleo 1C50 L29

The silicone sleeve - option

Cover below the silicone sleeve or without

Cover over the silicone sleeve

In this case, make sure you are taking the measurement of circuit with the sleeve on the socket.

Adapter

- 30
- 34

Other, e.g. Harmony P4:

Colour

- SkinTone
- Black
- Grey
- White

Skin 1..29
(charged extra)

Finish

- Matt varnish
(charged extra)

Tattoo (charged extra)

- Own
- Catalog number from the our website.
e.g.; 07-1001

Only for SkinTone, Red and White colour

Check list :

Photo of the foot from a distance of about 1 meter


camera lens position is very important

Camera lens placed perpendicular to the foot pyramid




Two photos of the prosthesis.

The camera lens is approximately 30 cm above the floor. There must not be any shoes, socks or other pieces of clothing included in the photo.



Scan of the prosthesis and the patient's leg *

in a standing position, without the shoes, socks and other items of clothing - *.stl; *.obj



- All our covers can be shortened from the bottom by a few cm.
- In the delivered condition, the lower part of the cover is flat.
- In the case of a parabolic foot shell, the lower part of the cover can be ground down to the shape of the foot.

