


## Circumferences of the prosthetic socket at heights of 50 mm



### Adapter

Above the knee

☐ 30

☐ 34

e.g.:

4WR95=3

### New Colours

21	41	61	81
23	43	63	83
25	45	65	85

Enter the number :



- matte finish as standard  
- not available for tattoos

### Colour

☐ Skin Tone1

☐ Skin Tone3

☐ Graphite

☐ White

☐ Grey

☐ Black

☐ Matt varnish (charged extra)

### Tattoo (charged extra)

Available only for Skin Tone3, White and Grey colour


☐ Catalog number  
e.g.: 08-1001

☐ Own

## Check list :

☐

**4 photos**, as shown in the images above. The camera must be placed perpendicular to the legs, in a low position, and both legs must be centered on the photo. Please remember that all prosthetic components must be completely visible, including the socket. Please remove shoes, socks and any object that covers the prosthesis.


☐

Scan of the prosthesis and the patient's leg



Format  
.stl  
.obj


☐

Scan of the entire prosthesis, without a patient.




Format  
.stl  
.obj



## NOTES :



- If the provided cover does not have a hole for the valve, and it is required, it should be made.
- After attaching the cover, the lower part should be approximately 5-6 cm below the top of the prosthetic knee.
- Detailed information on the installation of the cover can be found on the website and delivered together with the product.



NEW COLOURS - do not send.

