



The height of the cover in [mm]



We measure from the end of the knee (cover).

Prosthetic Socket circumference at the height in [mm]



0cm
-5cm
-10cm

0	
-5	
-10	

To take the measurement correctly, please mark the height of the cover in several places and place a ruler below the marking. The measurement should be made very carefully because moving the ruler even by a few millimetres up or down causes an incorrect reading of the circumference even by several centimetres.



Adapter
above the knee

- 30
- 34
- Lack

Colour

- SkinTone
- Black
- Grey
- White

Finish

- Matt varnish
(charged extra)

The matte version is not available for the black colour.

Tattoo (charged extra)

- Catalog number from our website.
e.g.; 08-1001

Only for SkinTone, Red and White colour

Check list :

Three photos of the prosthesis.

The camera lens should be at the height of the valve.



Scan of the prosthesis and the patient from the waist down

Before scanning, remove the prosthesis from the patient's stump so that the entire prosthetic socket is scanned. - *.stl; *.obj

- The hole for the valve or PIN should be made on your own, cutting and grinding to the desired size

