


Flex Neo AK

Flex Neo MAX AK

Type of amputation	Above Knee Transfemoral	Above Knee Transfemoral
Magnets	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Height	Any of the range of 32 - 53 cm	Any
Circumference	36; 38; 40; 42 cm	Any
Internal dimensions of the lower part	58 mm x 85 mm	According to the dimensions of the foot shells
Shape	Anatomical (left or right)	Reconstruction of calf based on a scan.
Cross-section of the lower part	Standard elliptical shape	Adapted to the shape of the shell of the foot
Colours		
Painless tattoo	<input type="checkbox"/>	Charged extra
Matte finish	Charged extra	Charged extra
A scan is required to order	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Cosmetic silicone 29 colours	<input type="checkbox"/>	<input checked="" type="checkbox"/>